

Heat Dissipates

- Extra Fine ●
- Fine ●
- Medium ●
- Coarse ●
- Super Coarse ●

Shape No.	368T	863T	836T	880T	850T	855T	850T	855T	856T	847T	847T	837KRT	850KRT
Size (1/10 mm)	023	014	012	010	014	016	018	018	018	016	018	012	018
Length (mm)	5.0	8.0	5.0	7.0	10.0	7.0	10.0	7.0	9.0	9.0	8.0	8.0	9.0
Order No.	285.5T	265.8T	515.5T	KS0T	771.10T	772.7T	772.10T	773.7T	773.9T	702.9T	703.8T	837KRT012	850KRT018

Occlusal Reduction / Finishing

Wheels

Needle

Mosquito

Modified Shoulder

Bulk Reducer

Depth Cutter 0.5 mm

Shape No.	379	368	368	368	811	909	909	889	858	859	852	392	837KR	850KR	6055	834
Size (1/10 mm)	023	016	018	023	033	035	043	008	014	014	023	016	012	018	018	021
Length (mm)	4.5	3.0	5.0	5.0	4.2	1.1	1.3	3.0	8.0	10.0	10.0	5.0	8.0	9.0	8.0	6.0
Order No.	275.5	292.3	280.5	285.5	230	860	863	201.3	1314.8	1314.10	209.10	392 016	837KR012	850KR018	ST8	834 021
Extra Fine 30µ				285.5XF				201.3XF		1314.10XF						
Fine 100µ	275.5F	292.3F	280.5	285.5F				201.3F	1314.8F	1314.10F	209.10F	392 016F				
Medium 130µ								201.3M	1314.8M	1314.10M		392 016M				834 021M
Coarse 145µ	275.5C	292.3C	280.5C	285.5C	230C	860C	863C		1314.8C	1314.10C	209.10C		837KR012C	850KR018C	ST8C	
Super Coarse 180µ	275.5SC	292.3SC	280.5SC	285.5SC	230SC	860SC	863SC								ST8SC	

Flame

Beveled Cylinder

Curettage

Round

Shape No.	879	862	863	863	862	863	863	878	879	878	878K	878K	878K	879K	801	801	801
Size (1/10 mm)	018	012	012	014	016	016	018	012	012	014	014	016	018	018	010	014	023
Length (mm)	9.0	8.0	10.0	8.0	8.0	10.0	9.0	8.0	10.0	8.0	9.0	9.0	8.0	11.0	1.0	1.4	2.3
Order No.	256.9	260.8	260.10	265.8	267.8	267.10	270.9	250.8	250.10	251.8	254SB	257SB	258SB	257SC	115	120	135
Extra Fine 30µ		260.8XF	260.10XF		267.8XF												120XF
Fine 100µ	256.9F	260.8F	260.10F	265.8F	267.8F	267.10F	270.9F	250.8F									135F
Medium 130µ		260.8M	260.10M		267.8M			250.8M	250.10M	251.8M	254SBM	257SBM	258SBM	257SCM	115M	120M	
Coarse 145µ	256.9C	260.8C	260.10C	265.8C	267.8C	267.10C	270.9C	250.8C	250.10C	251.8C	254SBC	257SBC	258SBC	257SCC	115C	120C	135C
Super Coarse 180µ	256.9SC			265.8SC				250.8SC									135SC

Inverted Cone

Double Cone

Flat End Cylinder

Pear

Round End Cylinder

Shape No.	805	805	805	813	806	835	835	835	835	836	837	830	830	838	880	880	881
Size (1/10 mm)	010	014	016	014	012	008	009	010	012	012	012	010	014	012	010	012	014
Length (mm)	1.0	2.0	2.0	2.0	2.5	4.0	3.0	4.0	4.0	5.0	7.0	2.2	3.0	4.0	7.0	7.0	8.0
Order No.	312	313	320.2	1220	1250	514.4	541	835-010	835-012	515.5	515.7	331	365.4	GW2	KS0	KS1	KS2
Medium 130µ	312M	313M		1220M	1250M			835-010M				331M	365.4M				
Coarse 145µ	312C	313C	320.2C			514.4C	541C	835-010C	835-012C	515.5C	515.7C	331C	365.4C	GW2C	KS0C	KS1C	KS2C
Super Coarse 180µ										515.5SC	515.7SC		365.4SC		KS0SC	KS1SC	KS2SC

Round End Taper

Flat End Taper

Shape No.	850	855	856	850	855	856	850	855	856	856	856	847	847	848	847	848	847	848
Size (1/10 mm)	012	012	014	014	016	016	018	018	018	021	023	012	014	016	016	018	018	023
Length (mm)	11/0.5	7/0.7	8/0.9	10/0.8	7/1.0	8/1.1	10/1.1	7/1.4	9/1.4	8/1.6	7/1.8	9/0.6	9/0.7	11/0.8	9/1.1	10/1.1	8/1.1	10/1.4
Order No.	770.11	771.7	771.8	771.10	772.7	772.8	772.10	773.7	773.9	774.8	775.7	700.9	701.9	701.11	702.9	702.10	703.8	703.10
Extra Fine 30µ	770.11XF																	
Fine 100µ	770.11F			771.10F		772.8F	772.10F	773.7F	773.9F			700.9F						
Medium 130µ	770.11M		771.8M	771.10M	772.7M		772.10M	773.7M	773.9M			701.9M		702.9M	702.10M	703.8M	703.10M	
Coarse 145µ	770.11C	771.7C	771.8C	771.10C	772.7C	772.8C	772.10C	773.7C	773.9C	774.8C	775.7C	700.9C	701.9C	701.11C	702.9C	702.10C	703.8C	703.10C
Super Coarse 180µ					772.7SC	772.8SC	772.10SC	773.7SC	773.9SC	774.8SC	775.7SC	700.9SC			702.9SC	702.10SC	703.8SC	703.10SC