

DENTAL LED HEADLIGHT

Features:

- Loupes mounted – Extra Light comes with a mount attachment for any Magnification Loupes.
- Daylight color light – allows the doctor to see tissues and teeth in the oral cavity in their true color, just like in natural day light.
- Narrow beam – allows the surgeon to illuminate the oral cavity, without dazzling the patient.
- Adjustable slope angle – beam’s angle can be adjusted up to 180 degrees to fit individual surgeons and tasks.
- Lightweight – Extra light weight just 15grams, applying minimal pressure on the forehead or the nose, when mounted on glasses.
- Cordless – Extra light is powered by one of two compact rechargeable power packs included, allowing the surgeon free of movement.

DENTAL LED HEADLIGHT



Item No. 08260A Snap on type complete set
08260B Clip on type complete set

Technical specifications

Headlight :

Luminous intensity: 130 lm (min)
 Brightness on 40 cm distance:
 40,000Lux
 Color temperature: 6000(±200) °K
 Illuminate beam radius on 40 cm
 distance: 40mm
 Dimensions (w/o headband): 30(L) x
 24(D) mm
 Weight (w/o headband): 15g
 Led lifetime: 50,000 hours

Item No. 08260A-4 Snap on type Headlight only
08260B-4 Clip on type only

Controller :

Input : 4.2V
 Battery capacity: 3,600mA
 Full battery charging time: 8 hours
 Average operating time: 8.5 hours
 Dimensions: 146(L) x 65(W) x 23(D)
 mm
 Weight: 142g

Item No. 08260A-3

Charger :

Type: AC/DC Adapter
 Input voltage: 100 – 240V
 Input frequency: 50 - 60 Hz
 Wall outlet plug: Euro / US
 Output Power: DC 4.2V / 1000mA

Item No. 08260A-2