

Szczecin 24 08 2011

Clinical Opinion

Re: Curing Light TC - 24 by Spring Health Products

Since using the Spring Health Curing Lights since 2009 I can state:

The Spring Health Curing Lights cure composites in a very short time, approximately half time of other lights. Most of the filling materials can be cured in 6 seconds.

The unbreakable, metal housing of the Curing Light is durable, easy to use and disinfect.

The Spring Health Curing Light is very easy to handle and very comfortable to use.


The Curing Light can be easily held in any position in the same manner as a turbine handpiece.

The Curing Light has two on/off switches, which increase the ergonomics when used in upper or lower arch.

The tested output power is over 1200 mW/cm². This makes it a very practical, ergonomic instrument.

This power makes it most efficient light I ever used.

Curing Light TC - 24 by Spring Health Products is extremely ergonomic thanks to its low weight and way of mounting into dental unit.

lek. stom. Paweł Andersz
ul. Okólna 43
71-742 Szczecin
tel. 501 652 932

lek. stom. Paweł Andersz